


TRPC6 Antibody (R33837)

Catalog No.	Formulation	Size
R33837-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Cow
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	7225
Applications	Western Blot : 0.03-0.1ug/ml ELISA (peptide) LOD : 1:128000
Limitations	This TRPC6 antibody is available for research use only.



Western blot testing of human placenta lysate with TRPC6 antibody at 0.03ug/ml.
Predicted molecular weight: ~106 kDa.

Description

Additional name(s) for this target protein: Transient receptor potential cation channel, subfamily C, member 6

Application Notes

Optimal dilution of the TRPC6 antibody should be determined by the researcher.

Immunogen

Amino acids EKLSMEPNQEETNR were used as the immunogen for this TRPC6 antibody.

Storage

Aliquot and store the TRPC6 antibody at -20oC.