


## Trpv2 Antibody (R35870)

Catalog No.	Formulation	Size
R35870-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Mouse, Rat
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	22368
Applications	Western Blot : 1-3ug/ml ELISA (peptide) LOD : 1:8000
Limitations	This Trpv2 antibody is available for research use only.



Western blot testing of rat brain lysate with Trpv2 antibody at 2ug/ml. Expected molecular weight: 85-95 kDa depending on glycosylation level.

## Description

Additional name(s) for this target protein: Transient receptor potential cation channel, subfamily V, member 2

## Application Notes

Optimal dilution of the Trpv2 antibody should be determined by the researcher.

## Immunogen

Amino acids SEEDHLPLQVLQSH were used as the immunogen for this Trpv2 antibody.

## Storage

Aliquot and store the Trpv2 antibody at -20oC.

## References (1)