

TRPV5 Antibody Biotin Conjugate (R35921BTN)

Catalog No.	Formulation	Size
R35921BTN	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Cow, Dog, Mouse, Pig, Rat
Format	Biotin Conjugate
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	56302
Applications	Western Blot : 1-3ug/ml ICC/IF : suitable ELISA (peptide) LOD : 1:16000
Limitations	This TRPV5 antibody is available for research use only.



Western blot testing of human pancreas lysate with biotinylated TRPV5 antibody at 1ug/ml. Predicted/observed molecular weight: ~85kDa.

Description

Additional name(s) for this target protein: Transient receptor potential cation channel, subfamily V, member 5

Application Notes

Optimal dilution of the TRPV5 antibody should be determined by the researcher.

Immunogen

Amino acids SHRGWEILRQNT were used as the immunogen for this TRPV5 antibody.

Storage

Aliquot and store the TRPV5 antibody at -20oC.