


TRX Antibody / Thioredoxin / TXN (RQ4619)


Catalog No.	Formulation	Size
RQ4619	0.5mg/ml if reconstituted with 0.2ml sterile DI water	100 ug

Bulk quote request


Availability	1-3 business days
Species Reactivity	Mouse, Rat
Format	Antigen affinity purified
Host	Rabbit
Clonality	Polyclonal (rabbit origin)
Isotype	Rabbit IgG
Purity	Antigen affinity purified
Buffer	Lyophilized from 1X PBS with 2% Trehalose and 0.025% sodium azide
UniProt	P10639
Localization	Nuclear, cytoplasmic
Applications	Western Blot : 0.5-1ug/ml Immunohistochemistry (FFPE) : 1-2ug/ml Immunocytochemistry (FFPE) : 1-2ug/ml Direct ELISA : 0.1-0.5ug/ml (mouse recombinant protein) Immunofluorescence : 2-4ug/ml
Limitations	This TRX antibody is available for research use only.




Western blot testing of mouse 1) testis, 2) kidney, 3) small intestine, 4) spleen, 5) stomach, 6) lung, 7) smooth muscle and 8) heart lysate with TXN antibody at 0.5ug/ml. Predicted molecular weight ~12 kDa.




Western blot testing of rat 1) testis, 2) kidney, 3) small intestine, 4) spleen, 5) stomach, 6) lung, 7) smooth muscle and 8) heart lysate with TRX antibody at 0.5ug/ml. Predicted molecular weight ~12 kDa.




IHC staining of FFPE mouse small intestine with TRX antibody at 1ug/ml. HIER: boil tissue sections in pH6, 10mM citrate buffer, for 10-20 min followed by cooling at RT for 20 min.



ICC staining of FFPE mouse LLC cells with TRX antibody at 1ug/ml. HIER: boil tissue sections in pH6, 10mM citrate buffer, for 10-20 min followed by cooling at RT for 20 min.



IHC staining of FFPE rat spleen with TRX antibody at 1ug/ml. HIER: boil tissue sections in pH6, 10mM citrate buffer, for 10-20 min followed by cooling at RT for 20 min.



Immunofluorescent staining of mouse NIH3T3 cells with TRX antibody (green, top left) and DAPI nuclear stain (blue, top right), and overlay (bottom).

Description

Thioredoxin is a class of small redox proteins known to be present in all organisms. It is mapped to 9q31.3. The protein encoded by this gene acts as a homodimer and is involved in many redox reactions. The encoded protein is active in the reversible S-nitrosylation of cysteines in certain proteins, which is part of the response to intracellular nitric oxide. This protein is found in the cytoplasm. Two transcript variants encoding different isoforms have been found for this gene.

Application Notes

Optimal dilution of the TRX antibody should be determined by the researcher.

Immunogen

Amino acids V2-A105 from the mouse protein were used as the immunogen for the TRX antibody.

Storage

After reconstitution, the TRX antibody can be stored for up to one month at 4°C. For long-term, aliquot and store at -20°C. Avoid repeated freezing and thawing.