


## TRX Antibody / TXN / Thioredoxin [clone AOGB-20] (RQ5327)

Catalog No.	Formulation	Size
RQ5327	Antibody in PBS with 0.02% sodium azide, 50% glycerol and 0.4-0.5mg/ml BSA	100 ul

[Bulk quote request](#)

Availability	1-2 weeks
Species Reactivity	Human
Format	Purified
Host	Rabbit
Clonality	Rabbit Monoclonal
Isotype	Rabbit IgG
Clone Name	AOGB-20
Purity	Affinity purified
UniProt	P10599
Applications	Western Blot : 1:1000-1:5000
Limitations	This TRX antibody is available for research use only.



Western blot testing of human HeLa cell lysate with TRX antibody. Predicted molecular weight ~12 kDa.

## Description

The protein encoded by the TXN gene acts as a homodimer and is involved in many redox reactions. The encoded protein is active in the reversible S-nitrosylation of cysteines in certain proteins, which is part of the response to intracellular nitric oxide. This protein is found in the cytoplasm. Two transcript variants encoding different isoforms have been found for this gene. [RefSeq]

## Application Notes

Optimal dilution of the TRX antibody should be determined by the researcher.

## Immunogen

A synthetic peptide specific to human TRX / Thioredoxin was used as the immunogen for the TRX antibody.

## Storage

Store the TRX antibody at -20oC.