


TSGA2 Antibody (R34230)

Catalog No.	Formulation	Size
R34230-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human, Mouse
Predicted Reactivity	Cow, Rat
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	89765
Applications	Western Blot : 0.1-0.3ug/ml ELISA (peptide) LOD : 1:4000
Limitations	This TSGA2 antibody is available for research use only.



Western blot testing of mouse testis lysate with TSGA2 antibody at 0.1ug/ml. Predicted molecular weight: ~34kDa.

Description

Additional name(s) for this target protein: Radial spoke head 1; RSPH1

Application Notes

Optimal dilution of the TSGA2 antibody should be determined by the researcher.

Immunogen

Amino acids HLNHRYQGKFLNKNP were used as the immunogen for this TSGA2 antibody.

Storage

Aliquot and store the TSGA2 antibody at -20oC.