


TSPAN14 Antibody (R35853)

Catalog No.	Formulation	Size
R35853-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Mouse
Predicted Reactivity	Human, Rat, Dog, Pig, Cow
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	81619
Applications	Western Blot : 1-3ug/ml ELISA (peptide) LOD : 1:128000
Limitations	This TSPAN14 antibody is available for research use only.



Western blot testing of mouse lung lysate with TSPAN14 antibody at 2ug/ml.
Predicted/observed molecular weight: 30 kDa.

Description

Additional name(s) for this target protein: Tetraspanin 14

Application Notes

Optimal dilution of the TSPAN14 antibody should be determined by the researcher.

Immunogen

Amino acids SDIEAVKAGHH were used as the immunogen for this TSPAN14 antibody.

Storage

Aliquot and store the TSPAN14 antibody at -20oC.