


## TSPAN5 Antibody (R35904)

Catalog No.	Formulation	Size
R35904-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human, Mouse, Rat
Predicted Reactivity	Cow, Dog, Pig
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	10098
Applications	Western Blot : 0.3-1ug/ml ELISA (peptide) LOD : 1:128000
Limitations	This TSPAN5 antibody is available for research use only.



Western blot testing of mouse brain lysate with TSPAN5 antibody at 0.3ug/ml. Predicted molecular weight: ~30kDa, routinely observed at 30~40kDa.

### Description

Additional name(s) for this target protein: Tetraspanin 5, NET-4

### Application Notes

Optimal dilution of the TSPAN5 antibody should be determined by the researcher.

## **Immunogen**

Amino acids SGKHYKGPEVSC were used as the immunogen for this TSPAN5 antibody.

## **Storage**

Aliquot and store the TSPAN5 antibody at -20oC.