


## TSPAN8 Antibody (R35436)

Catalog No.	Formulation	Size
R35436-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	7103
Applications	Western Blot : 0.01-0.03ug/ml ELISA (peptide) LOD : 1:16000
Limitations	This TSPAN8 antibody is available for research use only.



Western blot testing of human colon lysate with TSPAN8 antibody at 0.01ug/ml.  
Predicted molecular weight: ~26kDa.

## Description

Additional name(s) for this target protein: Tetraspanin 8, TM4SF3, transmembrane 4 superfamily member 3

## Application Notes

Optimal dilution of the TSPAN8 antibody should be determined by the researcher.

## Immunogen

Amino acids QRPCQSYNGKQVYKE were used as the immunogen for this TSPAN8 antibody.

## Storage

Aliquot and store the TSPAN8 antibody at -20oC.