


TSPAN8 Antibody (R35436)

Catalog No.	Formulation	Size
R35436-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	7103
Applications	Western Blot : 0.01-0.03ug/ml ELISA (peptide) LOD : 1:16000
Limitations	This TSPAN8 antibody is available for research use only.



Western blot testing of human colon lysate with TSPAN8 antibody at 0.01ug/ml.
Predicted molecular weight: ~26kDa.

Description

Additional name(s) for this target protein: Tetraspanin 8, TM4SF3, transmembrane 4 superfamily member 3

Application Notes

Optimal dilution of the TSPAN8 antibody should be determined by the researcher.

Immunogen

Amino acids QRPCQSYNGKQVYKE were used as the immunogen for this TSPAN8 antibody.

Storage

Aliquot and store the TSPAN8 antibody at -20oC.