

## TSPYL2 Antibody (R35852)

Catalog No.	Formulation	Size
R35852-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	64061
Applications	Western Blot : 0.3-1ug/ml ELISA (peptide) LOD : 1:64000
Limitations	This TSPYL2 antibody is available for research use only.



Western blot testing of Jurkat lysate with TSPYL2 antibody at 0.3ug/ml. Predicted molecular weight: ~79kDa.

## Description

Additional name(s) for this target protein: TSPY-like 2; Differentially Expressed Nucleolar TGF beta1 Target; DENTT

## Application Notes

Optimal dilution of the TSPYL2 antibody should be determined by the researcher.

## Immunogen

Amino acids SDIDETIHDIKISD were used as the immunogen for this TSPYL2 antibody.

## Storage

Aliquot and store the TSPYL2 antibody at -20oC.