


TSPYL2 Antibody (R35852)

Catalog No.	Formulation	Size
R35852-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	64061
Applications	Western Blot : 0.3-1ug/ml ELISA (peptide) LOD : 1:64000
Limitations	This TSPYL2 antibody is available for research use only.



Western blot testing of Jurkat lysate with TSPYL2 antibody at 0.3ug/ml. Predicted molecular weight: ~79kDa.

Description

Additional name(s) for this target protein: TSPY-like 2; Differentially Expressed Nucleolar TGF beta1 Target; DENTT

Application Notes

Optimal dilution of the TSPYL2 antibody should be determined by the researcher.

Immunogen

Amino acids SDIDETIHDIKISD were used as the immunogen for this TSPYL2 antibody.

Storage

Aliquot and store the TSPYL2 antibody at -20oC.