


TLL4 Antibody (R33637)

Catalog No.	Formulation	Size
R33637-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Cow, Mouse, Rat
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	9654
Applications	Immunohistochemistry (FFPE) : 5-10ug/ml ELISA (peptide) LOD : 1:32000
Limitations	This TLL4 antibody is available for research use only.



IHC testing of FFPE human bronchus with TLL4 antibody at 5ug/ml. HIER: steamed with pH6 citrate buffer, AP-staining.

Description

Additional name(s) for this target protein: Tubulin tyrosine ligase-like family, member 4

Application Notes

Optimal dilution of the TLL4 antibody should be determined by the researcher.

Immunogen

Amino acids DSQDPTKEIRFTE were used as the immunogen for this TTLL4 antibody.

Storage

Aliquot and store the TTLL4 antibody at -20oC.