


TTLL4 Antibody (R35779)

Catalog No.	Formulation	Size
R35779-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human, Mouse, Rat, Dog
Predicted Reactivity	Cow
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	9654
Applications	Western Blot : 1-3ug/ml ELISA (peptide) LOD : 1:16000
Limitations	This TTLL4 antibody is available for research use only.



Western blot testing of MDCK lysate with TTLL4 antibody at 1ug/ml. Predicted molecular weight: ~133kDa.

Description

Additional name(s) for this target protein: Tubulin tyrosine ligase-like family, member 4

Application Notes

Optimal dilution of the TTLL4 antibody should be determined by the researcher.

Immunogen

Amino acids RSSYRQKPYQQLE were used as the immunogen for this TTLL4 antibody.

Storage

Aliquot and store the TTLL4 antibody at -20oC.