


TTLL4 Antibody (R35779)

| Catalog No. | Formulation | Size |
|--------------|---|--------|
| R35779-100UG | 0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide | 100 ug |

[Bulk quote request](#)

| | |
|----------------------|--|
| Availability | 1-3 business days |
| Species Reactivity | Human, Mouse, Rat, Dog |
| Predicted Reactivity | Cow |
| Format | Antigen affinity purified |
| Host | Goat |
| Clonality | Polyclonal (goat origin) |
| Isotype | Goat Ig |
| Purity | Antigen affinity |
| Gene ID | 9654 |
| Applications | Western Blot : 1-3ug/ml ELISA (peptide) LOD : 1:16000 |
| Limitations | This TTLL4 antibody is available for research use only. |



Western blot testing of MDCK lysate with TTLL4 antibody at 1ug/ml. Predicted molecular weight: ~133kDa.

Description

Additional name(s) for this target protein: Tubulin tyrosine ligase-like family, member 4

Application Notes

Optimal dilution of the TTLL4 antibody should be determined by the researcher.

Immunogen

Amino acids RSSYRQKPYQQLE were used as the immunogen for this TTLL4 antibody.

Storage

Aliquot and store the TTLL4 antibody at -20oC.