


TUBB3 Antibody (R33581)

Catalog No.	Formulation	Size
R33581-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Cow, Mouse, Pig, Rat
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	10381
Applications	Western Blot : 1-3ug/ml ELISA (peptide) LOD : 1:1000
Limitations	This TUBB3 antibody is available for research use only.



Western blot testing of HepG2 lysate with TUBB3 antibody at 2ug/ml. Predicted molecular weight: ~50kDa.

Description

Additional name(s) for this target protein: Tubulin beta 3

Application Notes

Optimal dilution of the TUBB3 antibody should be determined by the researcher.

1. This TUBB3 antibody is specific for isoform 1 of the protein.

Immunogen

Amino acids DPSGNYVGSD were used as the immunogen for this TUBB3 antibody.

Storage

Aliquot and store the TUBB3 antibody at -20oC.