


TXNDC11 Antibody (R35347)

Catalog No.	Formulation	Size
R35347-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	51061
Applications	Western Blot : 1-3ug/ml ELISA (peptide) LOD : 1:16000
Limitations	This TXNDC11 antibody is available for research use only.



Western blot testing of human thymus lysate with TXNDC11 antibody at 1ug/ml.
Predicted molecular weight: ~108kDa.

Description

Additional name(s) for this target protein: Thioredoxin domain containing 11

Application Notes

Optimal dilution of the TXNDC11 antibody should be determined by the researcher.

Immunogen

Amino acids QLDYAEYVRRDSE were used as the immunogen for this TXNDC11 antibody.

Storage

Aliquot and store the TXNDC11 antibody at -20oC.