


## TXNDC4 Antibody (R33932)

Catalog No.	Formulation	Size
R33932-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Cow, Mouse, Rat
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	23071
Applications	Western Blot : 0.2-0.6ug/ml ELISA (peptide) LOD : 1:32000
Limitations	This TXNDC4 antibody is available for research use only.



Western blot testing of K562 lysate with TXNDC4 antibody at 0.2ug/ml. Predicted molecular weight: ~47kDa. Both observed bands are blocked by addition of immunizing peptide.

## Description

Additional name(s) for this target protein: Thioredoxin domain containing 4

## Application Notes

Optimal dilution of the TXNDC4 antibody should be determined by the researcher.

## Immunogen

Amino acids EQKDSDNYRVFER were used as the immunogen for this TXNDC4 antibody.

## Storage

Aliquot and store the TXNDC4 antibody at -20oC.