


## TXNDC5 Antibody (R35562)


Catalog No.	Formulation	Size
R35562-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

<b>Availability</b>	1-3 business days
<b>Species Reactivity</b>	Human, Mouse
<b>Predicted Reactivity</b>	Cow
<b>Format</b>	Antigen affinity purified
<b>Host</b>	Goat
<b>Clonality</b>	Polyclonal (goat origin)
<b>Isotype</b>	Goat Ig
<b>Purity</b>	Antigen affinity
<b>Gene ID</b>	81567
<b>Applications</b>	Western Blot : 0.1-0.3ug/ml ELISA (peptide) LOD : 1:16000
<b>Limitations</b>	This TXNDC5 antibody is available for research use only.



Western blot testing of human lymph node lysate with TXNDC5 antibody at 0.1ug/ml. Isoforms of ~48kDa and ~36kDa may be detected.



Western blot testing of mouse lymph node (A) and mouse intestine (B) lysates with TXNDC5 antibody at 0.1ug/ml. Predicted molecular weight (mouse): ~46kDa.

## Description

Additional name(s) for this target protein: Thioredoxin domain containing 5

## Application Notes

Optimal dilution of the TXNDC5 antibody should be determined by the researcher.

## Immunogen

Amino acids RDLESLREYVESQLQ were used as the immunogen for this TXNDC5 antibody. This sequence is common to human isoforms 1 and 2.

## Storage

Aliquot and store the TXNDC5 antibody at -20oC.