


TXNDC6 Antibody (R34003)

Catalog No.	Formulation	Size
R34003-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	347736
Applications	Western Blot : 2-5ug/ml ELISA (peptide) LOD : 1:8000
Limitations	This TXNDC6 antibody is available for research use only.



Western blot testing of human uterus lysate with TXNDC6 antibody at 2ug/ml. Predicted molecular weight: ~30kDa, routinely observed at ~38kDa.

Description

Additional name(s) for this target protein: Thioredoxin domain containing 6, TXL2

Application Notes

Optimal dilution of the TXNDC6 antibody should be determined by the researcher.

Immunogen

Amino acids ESSTQPRLKITDLD were used as the immunogen for this TXNDC6 antibody.

Storage

Aliquot and store the TXNDC6 antibody at -20oC.