


TZFP Antibody (R35721)

Catalog No.	Formulation	Size
R35721-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	27033
Applications	Western Blot : 2-5ug/ml ELISA (peptide) LOD : 1:8000
Limitations	This TZFP antibody is available for research use only.



Western blot testing of human bone marrow lysate with TZFP antibody at 2ug/ml.
Predicted molecular weight: ~53kDa.

Description

Additional name(s) for this target protein: Testis zinc finger protein, ZBTB32, FAZF

Application Notes

Optimal dilution of the TZFP antibody should be determined by the researcher.

Immunogen

Amino acids RPSTSPCCPSSSTT were used as the immunogen for this TZFP antibody.

Storage

Aliquot and store the TZFP antibody at -20oC.