


U2AF1L4 Antibody (R33731)

Catalog No.	Formulation	Size
R33731-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Dog, Mouse, Rat
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	199746
Applications	Western Blot : 0.1-0.3ug/ml ELISA (peptide) LOD : 1:32000
Limitations	This U2AF1L4 antibody is available for research use only.



Western blot testing of human cerebellum lysate with U2AF1L4 antibody at 0.1ug/ml.
Predicted molecular weight: ~21kDa.

Description

Additional name(s) for this target protein: U2 small nuclear RNA auxiliary factor 1-like 4

Application Notes

Optimal dilution of the U2AF1L4 antibody should be determined by the researcher.

Immunogen

Amino acids EDGERAVAELSNR were used as the immunogen for this U2AF1L4 antibody.

Storage

Aliquot and store the U2AF1L4 antibody at -20oC.