


UBA52 Antibody / 60S ribosomal protein L40 [clone 28U43] (RQ8346)

Catalog No.	Formulation	Size
RQ8346	Antibody in PBS with 0.02% sodium azide, 50% glycerol and 0.4-0.5mg/ml BSA	100 ul

Recombinant **RABBIT MONOCLONAL**

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Purified
Host	Rabbit
Clonality	Recombinant Rabbit Monoclonal
Isotype	Rabbit IgG
Clone Name	28U43
Purity	Affinity chromatography
UniProt	P62987
Applications	Western Blot : 1:500-1:2000
Limitations	This UBA52 antibody is available for research use only.



Western blot testing of human 293T cell lysate with UBA52 antibody. Predicted molecular weight ~15 kDa, commonly observed at 10-15 kDa.

Description

Ribosomal protein L40 is essential for translation of a subset of cellular transcripts, and especially for cap-dependent translation of vesicular stomatitis virus mRNAs. [UniProt]

Application Notes

Optimal dilution of the UBA52 antibody should be determined by the researcher.

Immunogen

A synthetic peptide derived from human UBA52 was used as the immunogen for the UBA52 antibody.

Storage

Store the UBA52 antibody at -20oC.