


UBASH3A Antibody (R33629)


Catalog No.	Formulation	Size
R33629-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	53347
Applications	Western Blot : 1-3ug/ml Immunohistochemistry (FFPE) : 3-5ug/ml ELISA (peptide) LOD : 1:8000
Limitations	This UBASH3A antibody is available for research use only.



Western blot testing of human spleen lysate with UBASH3A antibody at 1ug/ml.
Predicted molecular weight: ~70kDa.



IHC staining of FFPE human breast with UBASH3A antibody at 4ug/ml. HIER: steamed with pH6 citrate buffer, AP-staining.

Description

Additional name(s) for this target protein: Ubiquitin associated and SH3 domain containing A

Application Notes

Optimal dilution of the UBASH3A antibody should be determined by the researcher.

Immunogen

Amino acids DSGIRISSVFASP were used as the immunogen for this UBASH3A antibody.

Storage

Aliquot and store the UBASH3A antibody at -20oC.