


UBASH3A Antibody (R35407)

Catalog No.	Formulation	Size
R35407-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Dog
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	53347
Applications	Immunohistochemistry (FFPE) : 2.5-5ug/ml ELISA (peptide) LOD : 1:64000
Limitations	This UBASH3A antibody is available for research use only.



IHC testing of FFPE human spleen with UBASH3A antibody at 2.5ug/ml. HIER: steamed with pH6 citrate buffer, AP-staining.

Description

Additional name(s) for this target protein: Ubiquitin associated and SH3 domain containing A, STS2, Suppressor of T-cell receptor signaling 2, TULA

Application Notes

Optimal dilution of the UBASH3A antibody should be determined by the researcher.

Immunogen

Amino acids QQCSTPDGKYYRPD were used as the immunogen for this UBASH3A antibody.

Storage

Aliquot and store the UBASH3A antibody at -20oC.