


## UBC3B Antibody (R33159)

Catalog No.	Formulation	Size
R33159-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human, Rat
Predicted Reactivity	Mouse
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	54926
Applications	Immunohistochemistry (FFPE) : 2-4ug/ml Immunoelectron Microscopy : stains chromatoid bodies in rat spermatids ELISA (peptide) LOD : 1:4000
Limitations	This UBC3B antibody is available for research use only.



IHC testing of FFPE human testis with UBC3B antibody at 2ug/ml. HIER: steamed with pH6 citrate buffer, HRP-staining.

### Description

Additional name(s) for this target protein: Ubiquitin-conjugating enzyme UBC3B; UBE2R2

### Application Notes

Optimal dilution of the UBC3B antibody should be determined by the researcher.

## **Immunogen**

Amino acids AQQQMTSSQKALM were used as the immunogen for this UBC3B antibody.

## **Storage**

Aliquot and store the UBC3B antibody at -20oC.