


UBE2C Antibody (R35897)

Catalog No.	Formulation	Size
R35897-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Cow, Dog, Mouse, Rat
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	11065
Applications	Western Blot : 1-3ug/ml ELISA (peptide) LOD : 1:8000
Limitations	This UBE2C antibody is available for research use only.



Western blot testing of HEK293 lysate with UBE2C antibody at 1 ug/ml. Predicted molecular weight: ~19kDa.

Description

Additional name(s) for this target protein: Ubiquitin-conjugating enzyme E2C

Application Notes

Optimal dilution of the UBE2C antibody should be determined by the researcher.

1. This UBE2C antibody will detect isoforms 1/2/4.

Immunogen

Amino acids SGDKGISAFPESDN were used as the immunogen for this UBE2C antibody.

Storage

Aliquot and store the UBE2C antibody at -20oC.