


# **Ubiquitin Antibody (R30490)**


Catalog No.	Formulation	Size
R30490	0.5mg/ml if reconstituted with 0.2ml sterile DI water	100 ug

## **Bulk quote request**

Availability	1-3 business days
Species Reactivity	Human, Mouse, Rat
Format	Antigen affinity purified
Clonality	Polyclonal (rabbit origin)
Isotype	Rabbit IgG
Purity	Antigen affinity
Buffer	Lyophilized from 1X PBS with 2.5% BSA and 0.025% sodium azide/thimerosal
UniProt	P0CG47
Applications	Western Blot : 0.5-1ug/ml Immunohistochemistry (FFPE) : 0.5-1ug/ml
Limitations	This Ubiquitin antibody is available for research use only.



Western blot testing of Ubiquitin antibody and Lane 1: rat thymus; 2: human MCF-7; 3: MM231; 4: HeLa; 5: SMMC-7721 cell lysate.



IHC-P: Ubiquitin antibody testing of human lung cancer tissue

### **Description**

Ubiquitin(originally, ubiquitous immunopoietic polypeptide) is a small regulatory protein that has been found in almost all tissues(ubiquitously) of eukaryotic organisms. Its gene is mapped to 16p 13.3. The ubiquitin protein itself consists of 76 amino acids and has a molecular mass of about 8.5 kDa. Ubiquitin binds to proteins and labels them for destruction. It directs protein recycling and also direct proteins to other locations in the cell, where they control other protein and cell mechanisms.

### **Application Notes**

The stated application concentrations are suggested starting points. Titration of the Ubiquitin antibody may be required due to differences in protocols and secondary/substrate sensitivity.

### **Immunogen**

Amino acids 29-45 (KIQDKEGIPPDQQRLIF) were used as the immunogen for this Ubiquitin antibody (100% human, mouse and rat).

#### **Storage**

After reconstitution, the Ubiquitin antibody can be stored for up to one month at 4oC. For long-term, aliquot and store at -20oC. Avoid repeated freezing and thawing.