


## UBN1 Antibody (R33562)

Catalog No.	Formulation	Size
R33562-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Cow, Dog, Mouse, Pig, Rat, Zebrafish
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	29855
Applications	Western Blot : 1-3ug/ml ELISA (peptide) LOD : 1:128000
Limitations	This UBN1 antibody is available for research use only.



Western blot testing of MOLT4 lysate with UBN1 antibody at 1ug/ml. The expected ~120kDa band and the additional ~28kDa band are both blocked by the immunizing peptide.

## Description

Additional name(s) for this target protein: Ubinuclein 1

## Application Notes

Optimal dilution of the UBN1 antibody should be determined by the researcher.

## **Immunogen**

Amino acids DNSEAYDELVPASLT were used as the immunogen for this UBN1 antibody.

## **Storage**

Aliquot and store the UBN1 antibody at -20oC.