


## UBR4 Antibody (R34936)

Catalog No.	Formulation	Size
R34936-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Cow, Dog, Mouse, Rat
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	23352
Applications	Immunohistochemistry (FFPE) : 5-10ug/ml ELISA (peptide) LOD : 1:128000
Limitations	This UBR4 antibody is available for research use only.



IHC testing of FFPE human spleen with UBR4 antibody at 5ug/ml. HIER: steamed with pH6 citrate buffer, AP-staining.

### Description

Additional name(s) for this target protein: Ubiquitin protein ligase E3 component n recognin 4

### Application Notes

Optimal dilution of the UBR4 antibody should be determined by the researcher.

## Immunogen

Amino acids NTTRATSREEKNLQ were used as the immunogen for this UBR4 antibody.

## Storage

Aliquot and store the UBR4 antibody at -20oC.