


UBTD1 Antibody [clone 189CT14.5.4] (F40245)

Catalog No.	Formulation	Size
F40245-0.4ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.4 ml
F40245-0.08ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.08 ml

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Purified
Host	Mouse
Clonality	Monoclonal (mouse origin)
Isotype	Mouse IgG1, k
Clone Name	189CT14.5.4
Purity	Purified
UniProt	Q9HAC8
Applications	Western Blot : 1:100-1:500
Limitations	This UBTD1 antibody is available for research use only.



Description

Ubiquitin domain-containing protein 1

Application Notes

Titration of the UBTD1 antibody may be required due to differences in protocols and secondary/substrate sensitivity.

Immunogen

This UBTD1 antibody was produced from a mouse immunized with UBTD1 recombinant protein.

Storage

Aliquot the UBTD1 antibody and store frozen at -20oC or colder. Avoid repeated freeze-thaw cycles.