


## UCHL5 Antibody (F47733)

Catalog No.	Formulation	Size
F47733-0.4ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.4 ml
F47733-0.08ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.08 ml

[Bulk quote request](#)

<b>Availability</b>	1-3 business days
<b>Species Reactivity</b>	Human
<b>Format</b>	Antigen affinity purified
<b>Clonality</b>	Polyclonal (rabbit origin)
<b>Isotype</b>	Rabbit Ig
<b>Purity</b>	Antigen affinity
<b>UniProt</b>	Q9Y5K5
<b>Applications</b>	Western Blot : 1:1000
<b>Limitations</b>	This UCHL5 antibody is available for research use only.



## Description

Protease that specifically cleaves 'Lys-48'-linked polyubiquitin chains. Deubiquitinating enzyme associated with the 19S regulatory subunit of the 26S proteasome.

## Application Notes

Titration of the UCHL5 antibody may be required due to differences in protocols and secondary/substrate sensitivity.

## Immunogen

A portion of amino acids 136-164 from the human protein was used as the immunogen for this UCHL5 antibody.

## Storage

Aliquot the UCHL5 antibody and store frozen at -20oC or colder. Avoid repeated freeze-thaw cycles.