


UFM1 Antibody / Ubiquitin-fold modifier 1 [clone 28U42] (RQ8345)

Catalog No.	Formulation	Size
RQ8345	Antibody in PBS with 0.02% sodium azide, 50% glycerol and 0.4-0.5mg/ml BSA	100 ul

Recombinant **RABBIT MONOCLONAL**

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human, Mouse, Rat
Format	Purified
Host	Rabbit
Clonality	Recombinant Rabbit Monoclonal
Isotype	Rabbit IgG
Clone Name	28U42
Purity	Affinity chromatography
UniProt	P61960
Applications	Western Blot : 1:500-1:2000
Limitations	This UFM1 antibody is available for research use only.



Western blot testing of human HepG2 cell lysate with UFM1 antibody. Predicted molecular weight: 9-11 kDa (two isoforms).

Description

Ubiquitin-like modifier which can be covalently attached via an isopeptide bond to lysine residues of substrate proteins as a monomer or a lysine-linked polymer. [UniProt]

Application Notes

Optimal dilution of the UFM1 antibody should be determined by the researcher.

Immunogen

A synthetic peptide derived from human UFM1 was used as the immunogen for the UFM1 antibody.

Storage

Store the UFM1 antibody at -20oC.