


## UGT8 Antibody (R34802)

Catalog No.	Formulation	Size
R34802-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Cow, Dog, Rat
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	7368
Applications	Western Blot : 1-3ug/ml ELISA (peptide) LOD : 1:64000
Limitations	This UGT8 antibody is available for research use only.



Western blot testing of MCF7 overexpressing human UGT8 with UGT8 antibody at 1ug/ml [mock transfection in first lane].

## Description

Additional name(s) for this target protein: UDP glycosyltransferase 8

## Application Notes

Optimal dilution of the UGT8 antibody should be determined by the researcher.

## **Immunogen**

Amino acids NGKYKRNGHIK were used as the immunogen for this UGT8 antibody.

## **Storage**

Aliquot and store the UGT8 antibody at -20oC.