


UHRF1 Antibody (R34731)

Catalog No.	Formulation	Size
R34731-100UG	0.5~mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

Bulk quote request

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Dog, Mouse, Rat
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	29128
Applications	Western Blot : 1-3ug/ml ELISA (peptide) LOD : 1:2000
Limitations	This UHRF1 antibody is available for research use only.



Description

Additional name(s) for this target protein: Ubiquitin-like, containing PHD and RING finger domains, 1

Application Notes

Optimal dilution of the UHRF1 antibody should be determined by the researcher.

Immunogen

Amino acids LNQLFPGYGNGR were used as the immunogen for this UHRF1 antibody. **Storage** Aliquot and store the UHRF1 antibody at -20oC.