

UHRF1 Antibody (R34731)

Catalog No.	Formulation	Size
R34731-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Dog, Mouse, Rat
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	29128
Applications	Western Blot : 1-3ug/ml ELISA (peptide) LOD : 1:2000
Limitations	This UHRF1 antibody is available for research use only.



Western blot testing of nuclear Jurkat lysate with UHRF1 antibody at 2ug/ml.
Predicted/observed molecular weight: ~91/91-106kDa.

Description

Additional name(s) for this target protein: Ubiquitin-like, containing PHD and RING finger domains, 1

Application Notes

Optimal dilution of the UHRF1 antibody should be determined by the researcher.

Immunogen

Amino acids LNQLFPGYGNGR were used as the immunogen for this UHRF1 antibody.

Storage

Aliquot and store the UHRF1 antibody at -20oC.