


UHRF1BP1 Antibody (F53978)

Catalog No.	Formulation	Size
F53978-0.2ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.2 ml
F53978-0.05ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.05 ml

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Host	Rabbit
Clonality	Polyclonal (rabbit origin)
Isotype	Rabbit Ig
Purity	Antigen affinity
UniProt	Q6BDS2
Applications	Western Blot : 1:1000-2000
Limitations	This UHRF1BP1 antibody is available for research use only.



Western blot testing of human A431 cell lysate with UHRF1BP1 antibody at 1:2000.
Predicted molecular weight: 159 kDa, observed here at ~80 kDa.

Description

Ubiquitin-like containing PHD and RING finger domains 1-binding protein 1 may act as a negative regulator of cell growth.

Application Notes

Titration of the UHRF1BP1 antibody may be required due to differences in protocols and secondary/substrate sensitivity.

Immunogen

A portion of amino acids 1096-1129 from the human protein was used as the immunogen for the UHRF1BP1 antibody.

Storage

Aliquot the UHRF1BP1 antibody and store frozen at -20oC or colder. Avoid repeated freeze-thaw cycles.