


ULK3 Antibody (R34533)

Catalog No.	Formulation	Size
R34533-100UG	0.5~mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

Bulk quote request

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Dog, Pig
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	25989
Applications	Western Blot : 1-3ug/ml ELISA (peptide) LOD : 1:32000
Limitations	This ULK3 antibody is available for research use only.



Description

Additional name(s) for this target protein: Unc-51-like kinase 3

Application Notes

Optimal dilution of the ULK3 antibody should be determined by the researcher.

Immunogen

Amino acids KMRESRWEADTLDK were used as the immunogen for this ULK3 antibody. **Storage** Aliquot and store the ULK3 antibody at -20oC.