


ULK3 Antibody (R34533)

Catalog No.	Formulation	Size
R34533-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Dog, Pig
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	25989
Applications	Western Blot : 1-3ug/ml ELISA (peptide) LOD : 1:32000
Limitations	This ULK3 antibody is available for research use only.



Western blot testing of human hippocampus lysate with ULK3 antibody at 1ug/ml. Predicted molecular weight: ~53kDa. Both observed bands are blocked by addition of immunizing peptide.

Description

Additional name(s) for this target protein: Unc-51-like kinase 3

Application Notes

Optimal dilution of the ULK3 antibody should be determined by the researcher.

Immunogen

Amino acids KMRESRWEADTLDK were used as the immunogen for this ULK3 antibody.

Storage

Aliquot and store the ULK3 antibody at -20oC.