


UNC5CL Antibody (R34456)

Catalog No.	Formulation	Size
R34456-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	222643
Applications	Western Blot : 0.3-1ug/ml ELISA (peptide) LOD : 1:8000
Limitations	This UNC5CL antibody is available for research use only.



Western blot testing of human kidney lysate with UNC5CL antibody at 0.3ug/ml.
Predicted molecular weight: ~58kDa.

Description

Additional name(s) for this target protein: Unc-5 homolog C (C. elegans)-like

Application Notes

Optimal dilution of the UNC5CL antibody should be determined by the researcher.

Immunogen

Amino acids KHCAEQPSHARTY were used as the immunogen for this UNC5CL antibody.

Storage

Aliquot and store the UNC5CL antibody at -20oC.