


uPA Antibody / Urokinase / PLAU [clone AABC-16] (RQ4847)

Catalog No.	Formulation	Size
RQ4847	Antibody in PBS with 0.02% sodium azide, 50% glycerol and 0.4-0.5mg/ml BSA	100 ul

[Bulk quote request](#)

Availability	1-2 weeks
Species Reactivity	Human
Format	Purified
Host	Rabbit
Clonality	Rabbit Monoclonal
Isotype	Rabbit IgG
Clone Name	AABC-16
Purity	Affinity purified
UniProt	P00749
Applications	Western Blot : 1:500-1:1000
Limitations	This uPA antibody is available for research use only.



Western blot testing of human LNCaP cell lysate with uPA antibody. Predicted molecular weight ~48 kDa.

Description

'Urokinase-type plasminogen activator' specifically cleaves the zymogen plasminogen to form the active enzyme plasmin. [UniProt]

Application Notes

Optimal dilution of the uPA antibody should be determined by the researcher.

Immunogen

A synthetic peptide specific to human Urokinase was used as the immunogen for the uPA antibody.

Storage

Store the uPA antibody at -20oC.