


UPAR Antibody (R35785)

Catalog No.	Formulation	Size
R35785-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	5329
Applications	Western Blot : 1-3ug/ml ELISA (peptide) LOD : 1:128000
Limitations	This UPAR antibody is available for research use only.



Western blot testing of U251 lysate with UPAR antibody at 1ug/ml. Predicted molecular weight: ~37/60kDa (unmodified/glycosylated).

Description

Additional name(s) for this target protein: U-plasminogen activator receptor; PLAUR; CD87

Application Notes

Optimal dilution of the UPAR antibody should be determined by the researcher.

Immunogen

Amino acids RTGLKITSLTE were used as the immunogen for this UPAR antibody.

Storage

Aliquot and store the UPAR antibody at -20oC.