


Urocortin 3 Antibody (R34575)

Catalog No.	Formulation	Size
R34575-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	114131
Applications	Western Blot : 1-3ug/ml ELISA (peptide) LOD : 1:16000
Limitations	This Urocortin 3 antibody is available for research use only.



Western blot testing of HEK293 lysate over expressing human UCN3-FLAG with Urocortin 3 antibody [1ug/ml] in Lane A and probed with anti-FLAG (1/3000) in lane C. Mock-transfected HEK293 probed with Urocortin 3 antibody [1ug/ml] in Lane B. Predicted molecular weight: ~18kDa.

Description

Additional name(s) for this target protein: SCP, SPC, UCNIII; UCN3

Application Notes

Optimal dilution of the Urocortin 3 antibody should be determined by the researcher.

Immunogen

Amino acids KPRQDTAKSPHRTK were used as the immunogen for this Urocortin 3 antibody.

Storage

Aliquot and store the Urocortin 3 antibody at -20oC.