


UROD Antibody (R34358)

Catalog No.	Formulation	Size
R34358-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	7389
Applications	Western Blot : 0.1-0.3ug/ml ELISA (peptide) LOD : 1:128000
Limitations	This UROD antibody is available for research use only.



Western blot testing of K562 lysate with UROD antibody at 0.1ug/ml. The expected ~40kDa band and the additional ~28kDa band are both blocked by the immunizing peptide.

Description

Additional name(s) for this target protein: Uroporphyrinogen decarboxylase

Application Notes

Optimal dilution of the UROD antibody should be determined by the researcher.

Immunogen

Amino acids KARECVGKTVTLQGN were used as the immunogen for this UROD antibody.

Storage

Aliquot and store the UROD antibody at -20oC.