


Uromodulin Antibody (R33606)

Catalog No.	Formulation	Size
R33606-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	7369
Applications	Western Blot : 0.5-1.5ug/ml ELISA (peptide) LOD : 1:128000
Limitations	This Uromodulin antibody is available for research use only.



Western blot testing of HeLa lysate with Uromodulin antibody at 0.5ug/ml. Expected molecular weight: 70-105 kDa depending on level of glycosylation.

Description

Additional name(s) for this target protein: UMOD

Application Notes

Optimal dilution of the Uromodulin antibody should be determined by the researcher.

Immunogen

Amino acids DRCPHTRDSTIQ were used as the immunogen for this Uromodulin antibody.

Storage

Aliquot and store the Uromodulin antibody at -20oC.

References (1)