


# **Uromodulin Antibody (R35032)**

Catalog No.	Formulation	Size
R35032-100UG	0.5~mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

### **Bulk quote request**

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	7369
Applications	Western Blot : 0.03-0.1ug/ml ELISA (peptide) LOD : 1:128000
Limitations	This Uromodulin antibody is available for research use only.



### **Description**

Additional name(s) for this target protein: ADMCKD2, FJHN, HNFJ, HNFJ1, Tamm-Horsfall glycoprotein; UMOD

#### **Application Notes**

Optimal dilution of the Uromodulin antibody should be determined by the researcher.

#### **Immunogen**

Amino acids PITRKGVQATVS were used as the immunogen for this Uromodulin antibody.

## **Storage**

Aliquot and store the Uromodulin antibody at -20oC.

References (1)