


USH1C Antibody (R33613)

Catalog No.	Formulation	Size
R33613-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Cow, Mouse, Rat
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	10083
Applications	Western Blot : 0.1-0.3ug/ml ELISA (peptide) LOD : 1:128000
Limitations	This USH1C antibody is available for research use only.



Western blot testing of HEK293 lysate with USH1C antibody at 0.1ug/ml. Predicted molecular weight: ~62kDa, observed here at ~75kDa.

Description

Additional name(s) for this target protein: Usher syndrome 1C; Harmonin

Application Notes

Optimal dilution of the USH1C antibody should be determined by the researcher.

Immunogen

Amino acids DRKVAREFRHKVD were used as the immunogen for this USH1C antibody.

Storage

Aliquot and store the USH1C antibody at -20oC.