


## USP11 Antibody (R33347)

Catalog No.	Formulation	Size
R33347-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

<b>Availability</b>	1-3 business days
<b>Species Reactivity</b>	Human
<b>Format</b>	Antigen affinity purified
<b>Host</b>	Goat
<b>Clonality</b>	Polyclonal (goat origin)
<b>Isotype</b>	Goat Ig
<b>Purity</b>	Antigen affinity
<b>Gene ID</b>	8237
<b>Applications</b>	Western Blot : 1-3ug/ml ELISA (peptide) LOD : 1:8000
<b>Limitations</b>	This USP11 antibody is available for research use only.



Western blot testing of Jurkat lysate with USP11 antibody at 1ug/ml. Predicted molecular weight: ~110kDa, routinely observed at ~120kDa. (1)

### Description

Additional name(s) for this target protein: Ubiquitin specific protease 11

### Application Notes

Optimal dilution of the USP11 antibody should be determined by the researcher.

1. Possible degradation products of 55, 50 and 44 kDa may also be seen. (Ref. 1)

## Immunogen

Amino acids CSSPPSSEFMDVN were used as the immunogen for this USP11 antibody.

## Storage

Aliquot and store the USP11 antibody at -20oC.

## References (1)