

USP11 Antibody (R33347)

Catalog No.	Formulation	Size
R33347-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	8237
Applications	Western Blot : 1-3ug/ml ELISA (peptide) LOD : 1:8000
Limitations	This USP11 antibody is available for research use only.



Western blot testing of Jurkat lysate with USP11 antibody at 1ug/ml. Predicted molecular weight: ~110kDa, routinely observed at ~120kDa. (1)

Description

Additional name(s) for this target protein: Ubiquitin specific protease 11

Application Notes

Optimal dilution of the USP11 antibody should be determined by the researcher.

1. Possible degradation products of 55, 50 and 44 kDa may also be seen. (Ref. 1)

Immunogen

Amino acids CSSPPSSEFMDVN were used as the immunogen for this USP11 antibody.

Storage

Aliquot and store the USP11 antibody at -20oC.

References (1)