


USP15 Antibody (R34776)

Catalog No.	Formulation	Size
R34776-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Cow, Dog, Mouse, Pig, Rat
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	9958
Applications	Immunohistochemistry (FFPE) : 2.5-5ug/ml ELISA (peptide) LOD : 1:2000
Limitations	This USP15 antibody is available for research use only.



IHC testing of FFPE human liver with USP15 antibody at 2ug/ml. HIER: steamed with pH6 citrate buffer, HRP-staining.

Description

Additional name(s) for this target protein: Ubiquitin specific protease 15

Application Notes

Optimal dilution of the USP15 antibody should be determined by the researcher.

Immunogen

Amino acids NDIENENCMHTN were used as the immunogen for this USP15 antibody.

Storage

Aliquot and store the USP15 antibody at -20oC.